

FOR ANIMAL TREATMENT ONLY  
RESTRICTED VETERINARY  
MEDICINE

 Observe™  
avian

**TUBERCULIN**  
25,000 IU per mL

Registered to:ASUREQuality Limited,  
Level 1, 7a Pacific Rise, Mt Wellington, 1060  
New Zealand  
Tel: +64 9 573 8000  
Fax: +64 9 573 8118

Registered pursuant to ACVM Act 1997  
No. A10883.

See [www.foodsafety.govt.nz](http://www.foodsafety.govt.nz) for  
registration conditions

 ASUREQuality  
Kaitiaki Kai

STORE AT 2-8 DEGREES CELSIUS  
DO NOT FREEZE  
PROTECT FROM LIGHT  
STORE OUT OF THE REACH OF CHILDREN

FOR ANIMAL TREATMENT ONLY  
RESTRICTED VETERINARY MEDICINE

 Observe™  
avian

**TUBERCULIN**  
25,000 IU per mL

**Active Ingredients:** Purified Protein Derivative  
of *Mycobacterium avium*, strain D4ER,  
25,000 IU/mL.

**Indications:** For use as an aid in the diagnosis  
of mycobacterial infection in cattle and deer.

**Directions for use:** For use under the authority  
or prescription of a veterinary surgeon or by  
persons approved under Section 103 of the  
Biosecurity Act 1993 when applying cattle and  
deer tuberculin tests. Shake well before use.

#### Test Application

The following instructions are guidelines only:  
Dose: 0.1mL intra-dermally. It is recommended  
syringes are calibrated to 0.1 mL, and the  
needle of the syringe is 22-26 gauge and 3-4  
mm long. Syringes should be kept clean and  
sterilized. Contamination of the syringe with  
disinfectants or alcohol may interfere with  
testing. When filling syringes, care must be  
taken to prevent contamination of vials.

#### Preferred site of injection:

**Cattle:** the mid cervical site as a comparative  
cervical test (CCT);

**Deer:** the mid cervical site as a CCT.

When using the cervical sites, care is required  
to evenly clip hair close to the skin surface  
(2mm or less mean length) prior to injecting.  
The recommended size of the clipped area is  
10 x 10 cm for each injection site.

When applying a CCT use Observe™ Bovine  
Tuberculin 30,000 IU per mL in conjunction with  
Observe™ Avian Tuberculin 25,000 IU per mL.

#### Test Interpretation:

It is recommended the test be read 72 hours  
after injection.

Skin test reactions are to be classified in  
accordance with the requirements of OSPRI and

Biosecurity (National Bovine Tuberculin Pest  
Management Plan) Order 1998.

**Retesting:** To avoid desensitization, a retest  
using bovine and/or avian PPD tuberculin  
should not be applied until 60 days after any  
previous tuberculin test.

#### Product Information

**Withholding period:** Nil

**Storage:** Refrigerate (2°C – 8°C). Do not freeze.

Keep the vials in the closed box in order to  
protect from light.

**Transportation:** May be transported at 2°C –  
37°C for a period not longer than 14 days.  
Do not freeze.

**Use:** Do not use after the expiry date. Use  
immediately after opening. Discard unused  
contents. Do not mix with other vaccines or  
immunological products.

**Method of Disposal:** Unused, part-used and  
empty vials should be destroyed according to  
local regulations.

**Safety Precautions:** In the case of accidental  
self-injection, an area of intense local irritation  
may develop for up to 96 hours later. When  
irritation occurs, you are advised to consult your  
medical practitioner.

Store out of reach of children.

**Manufactured by:** Prionics Lelystad B.V.,  
Platinastraat 33, 8211 AR Lelystad,  
The Netherlands

Registered to: ASUREQuality Limited,  
Level 1, 7a Pacific Rise, Mt Wellington, 1060,  
New Zealand  
Tel: +64 9 573 8000  
Fax: +64 9 573 8118  
Registered pursuant to ACVM Act 1997  
No. A10883. See [www.foodsafety.govt.nz](http://www.foodsafety.govt.nz) for  
registration conditions.